

Product Highlight

Company Name

U.P. MARKETING GENERAL SUPPLY CO., LTD.

Product Highlight 1



Description Product
Highlight 1

IKA Bomb Calorimeter C6000 Global

C 6000 global standards oxygen bomb calorimeter combines modern technology, variability and automation in one instrument.

It operates according to all bomb calorimeter standards, such as e.g. the DIN, ISO, ASTM, GOST and GB.

The operator can choose between three different starting temperatures 22°C, 25°C, 30°C in each measuring mode.

Product Highlight 2



Description Product Highlight 2

Memrapure Aracus Amino acid analyzer

ARACUS advanced with two-pump configuration for independent operation of two fluidic lines (eluent side and reactive side)

96 samples Continuous analysis in a cooled autosampler

Post column derivatization with ninhydrin method

Reproducibility of retention time & area < 0.3 % RSD ca. 1 % RSD

Detection limit < 2 nmol/mL

Product Highlight 3



Description Product Highlight 3

DAIHAN ULT Freezer Fre800-B90

-90°C to -65°C SMART Multiuse Ultra Low Temperature Freezer with DuoFreez System
Remote Controlled by Mobile-Phone with WiRe™ App Service
High-quality Insulation Panel and Inner Doors ensure Defrosting around Outer Surface
Automatic Vacuum Breaker : Easy to re-open the Door
Capacity 796 Liter
CFC Free Mixed Gas Refrigerant

Product Highlight 4



Description Product Highlight 4

Aquinity² P35 Ultra Pure Water System

Produces Pure Water and Ultra Pure Water

Compact housing

Big touchscreen

Remote service module

Menus and wizards easy to use and mistake avoiding

Fixed and flexible dispenser

Easy access for cartridge exchange

Product Highlight 5



Description Product Highlight 5

Smart Automatic Door Analytical Balance model AP225W-AD

Capacity 220g

Display OEL display (dot matrix)

Pan Size \varnothing 91 mm

Sensitivity Stability Against Temperature Range ± 2 ppm/ $^{\circ}$ C (10 to 30 $^{\circ}$ C)

Ionizer Included

Adjustable Internal Windbreak Plate Included

Touchless Sensors Included

Automatic Doors Included

Clock-CAL Included

Sample (Concentration) Preparation Included

Minimum display 0.01 mg

mol Conversion Function Included

HPLC Buffer Solution Preparation Included

Recipe Compounding Included